Immuli<sup>TM</sup>

# Celiac tTG

rHuman Tissue Transglutaminase Antibody ELISA

**IVD** For *in vitro* diagnostic use

## PRODUCT INSERT

<b>REF</b>	5144A	Celiac tTG IgA ELISA	96 Determinations
<b>REF</b>	5144G	Celiac tTG IgG ELISA	96 Determinations

## INTENDED USE

Enzyme linked immunoassays (ELISA) for the qualitative and semi-quantitative detection of anti-human Tissue Transglutaminase IgA or IgG antibodies in human serum to aid in the diagnosis of gluten sensitive enteropathy / celiac disease (CD) in conjunction with other laboratory and clinical findings.

## SUMMARY AND EXPLANATION

Celiac Disease (CD) is an autoimmune gastrointestinal disorder that may occur in genetically susceptible individuals triggered by the ingestion of gluten-containing grains such as wheat, barley and rye. The classical symptoms of CD include diarrhea, weight loss and malnutrition. Only a small percentage of patients with CD presents with classical symptoms. Consequently, the clinical spectrum of CD has grown much broader than in the past to include patients that do not present with classical symptoms. It is not uncommon for the initial symptoms to be non-gastrointestinal or for gastrointestinal symptoms, if present, to be mild or intermittent. The need to examine a wider range of clinical presentation has led to greater numbers of individuals diagnosed with CD later in life than ever before. Adults may present with iron deficiency, macrocytic anemia and hypocalcaemia.

Studies have found the prevalence of CD to be highly variable. If clinical criteria alone are used in determining prevalence, the incidence of CD is much lower as compared with incidence established by serological methods.<sup>1,2</sup> Using serological methods, recent studies indicate the incidence of CD in the general population to be between one in 100 and one in 500.

Failure to diagnose CD early on may predispose an individual to long-term complications such as splenic atrophy and intestinal lymphoma.<sup>3,4</sup> A gluten-free diet (GFD) normalizes the mucosa and helps reduce the malignant potential.

Histological examination of the small intestinal biopsy remains the gold standard for diagnosing CD, but it has its own limitations. These include some patients with latent or even active CD that may have normal histopathology.<sup>5</sup>

The revised European Society of Paediatric Gastroenterology and Nutrition (ESPGHAN) criteria include only a single biopsy with clear-cut remission of clinical symptoms on GFD.<sup>6</sup> Positive serology at the time of diagnosis with

disappearance on GFD contributes to the diagnosis. The various serological tests employed in the work-up of patients suspected to have CD include gliadin (AGA), endomysial (EMA), and tissue transglutaminase (tTG) antibody tests. Antibodies to gliadin and tTG are detected by ELISA, whereas EMA are detected by indirect immunofluorescence. EMA are very specific indicators of CD. However, the EMA test is an immunohistochemical method that requires experience in reading immunofluorescence reactions.<sup>7</sup>

Since identification of tTG as the endomysial antigen, ELISA methods have been described for detecting antibodies in the sera of patients with CD. The advantage of the anti-tTG antibody assay is that it is automatable and less subjective than EMA. For this reason, many laboratories have opted to use the tTG antibody method as the screening method. Immulisa™ Celiac tTG is a unique immunoassay utilizing special chemistry of detecting antibodies to tTG with great degree of specificity and sensitivity and minimizes non-specific antibody detection in other disease controls as has been reported in some other tTG immunoassays.<sup>8-11</sup>

## PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with recombinant human Tissue Transglutaminase antigen followed by a blocking step to reduce non-specific binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the tTG antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgA or IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml) and reported as positive or negative.

## REAGENTS

### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

## Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.<sup>12</sup>

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

## Materials provided

ImmuliS<sup>TM</sup> Celiac tTG IgA ELISA

**REF** 5144A

ImmuliS<sup>TM</sup> Celiac tTG IgG ELISA

**REF** 5144G

Kits contains sufficient reagents to perform 96 determinations.

12 x 8	<b>MICROPLATE C tTG</b>	Microplate with individual breakaway microwells. Coated with recombinant human tTG. Ready for use.
1 x 1.75 ml	<b>CONTROL + C tTG-A</b>	Ready to use <b>Positive Control</b> (red cap) for <b>REF</b> 5144A. Contains human serum positive for tTG IgA antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	<b>CONTROL + C tTG-G</b>	Ready to use <b>Positive Control</b> (red cap) for <b>REF</b> 5144G. Contains human serum positive for tTG IgG antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	<b>CONTROL -</b>	Ready to use <b>Negative Control</b> (white cap). Contains human serum.
5 x 1.75 ml	<b>CALIBRATOR A C tTG-A</b> <b>CALIBRATOR B C tTG-A</b> <b>CALIBRATOR C C tTG-A</b> <b>CALIBRATOR D C tTG-A</b> <b>CALIBRATOR E C tTG-A</b>	Ready to use set of 5 <b>Calibrators</b> for <b>REF</b> 5144A. Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing tTG IgA antibodies. Concentrations in EU/ml are printed on the labels.
5 x 1.75 ml	<b>CALIBRATOR A C tTG-G</b> <b>CALIBRATOR B C tTG-G</b> <b>CALIBRATOR C C tTG-G</b> <b>CALIBRATOR D C tTG-G</b> <b>CALIBRATOR E C tTG-G</b>	Ready to use set of 5 <b>Calibrators</b> for <b>REF</b> 5144G. Calibrator A (green cap) 320 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing tTG IgG antibodies. Concentrations in EU/ml are printed on the labels.
1 x 15 ml	<b>IgA-CONJ HRP</b>	HRP goat anti-human <b>IgA Conjugate</b> for <b>REF</b> 5144A. Ready for use. Color coded pink.
1 x 15 ml	<b>IgG-CONJ HRP</b>	HRP goat anti-human <b>IgG Conjugate</b> for <b>REF</b> 5144G. Ready for use. Color coded pink.
1 x 60 ml	<b>DIL</b>	Serum Diluent. Ready for use. Color coded purple.

1 x 15 ml      SUBSTRATE|TMB

TMB enzyme substrate. Ready for use. **Protect from light.**

1 x 15 ml      STOP|H<sub>2</sub>SO<sub>4</sub>

Stop Solution. Ready for use.

2 x vials      BUF|WASH

Powder Wash Buffer. **Reconstitute to one liter each.**

1 x

Protocol Sheets

### Optional Components

1 x 60ml      BUF|WASH

Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

### Symbols used on labels

Lot number

Catalog number

In vitro diagnostic use

Use by

Storage temperature

Read instructions before use

Number of tests

Manufacturer

### Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

# PROCEDURE

## Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

## Test Method

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.
- or**
- For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative			
A	Blank	S5	Etc.
B	-	S6	
C	+	S7	
D	Cal D	S8	
E	S1	S9	
F	S2	S10	
G	S3	S11	
H	S4	S12	
	1	2	3

Semi-Quantitative			
A	Blank	S1	Etc.
B	-	S2	
C	+	S3	
D	Cal A	S4	
E	Cal B	S5	
F	Cal C	S6	
G	Cal D	S7	
H	Cal E	S8	
	1	2	3

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500µl** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

## Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be  $<10$  EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

## RESULTS

### Calculations

The concentrations of the patient samples can be determined by either of two methods:

## 1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

It is recommended that qualitative results be reported as “positive” or “negative.” Sample results greater than or equal to Calibrator D are considered positive.

## 2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with EU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods, such as assays for detection of EMA and/or Gliadin antibodies.

### **Interpretation**

The following serves only as a guide in the interpretation of the laboratory results. These values were determined by testing 114 normal blood donors and non-celiac disease control specimens. The mean of the normal subjects plus 2 SD was established as the assay cutoff and assigned an arbitrary value of 20 EU/ml. Each laboratory must determine its own normal values.

<u>anti-tTG Ab value</u>	<u>Interpretation</u>
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

### **Calibrator**

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

## **LIMITATIONS OF PROCEDURE**

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2-8°C for no

longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. Test results serve as an aid in the diagnosis and should be considered in conjunction with other laboratory and clinical findings.

**EXPECTED VALUES**

Test results in a normal population are usually expected to be negative. However, as the incidence of CD in the normal population is about 1%, some apparently healthy, asymptomatic individuals may test positive for tTG antibodies.

The incidence and levels of anti-tTG antibodies are dependent upon the diet status. The levels of these antibodies decrease and eventually will become negative in patients with CD who are on a gluten-free diet. Similarly, the levels of these antibodies will increase and may become positive when patients with CD who were on a gluten-free diet ingest a gluten-containing diet.<sup>13-15</sup> Patients who have CD but are IgA deficient will be positive for IgG antibodies to tTG. In such cases studies can be performed to confirm that the patient is IgA deficient.

**PERFORMANCE CHARACTERISTICS**

The utility of the Immulisa™ Celiac tTG ELISAs was evaluated by testing well-characterized EMA positive serum specimens from suspected CD subjects alongside disease controls and “normal” human sera. These specimens were also tested on commercially available ELISA and immunofluorescence test kits. These results are summarized below.

A. Immulisa™ Celiac tTG ELISAs vs. other tTG immunoassay:

		<b>Other tTG IgA ELISA</b>		
		Positive	Negative	Total
<b>IMMCO</b>	Positive	74	19	93
<b>CELIAC tTG</b>	Negative	2	90	92
<b>IgA ELISA</b>	Total	76	109	185

Positive Percent Agreement: 97.4% (95% CI 90.0% to 99.5%)

Negative Percent Agreement: 82.6% (95% CI 73.9% to 88.9%)

Overall Percent Agreement: 88.6% (95% CI 83.0% to 92.7%)

EMA Positive Celiac Subjects: 93

Disease Controls: 31

Healthy Normal Subjects: 61

### Other tTG IgG ELISA

		Positive	Negative	Total
<b>IMMCO</b>	Positive	74	21	95
<b>CELIAC tTG</b>	Negative	31	184	215
<b>IgG ELISA</b>	Total	105	205	310

Positive Percent Agreement: 70.5% (95% CI 60.7% to 78.8%)

Negative Percent Agreement: 89.8% (95% CI 84.6% to 93.4%)

Overall Percent Agreement: 83.2% (95% CI 78.5% to 87.1%)

EMA Positive Celiac Subjects: 166

Disease Controls: 53

Healthy Normal Subjects: 91

B. Immulisa™ Celiac tTG ELISAs vs. EMA: Results obtained with the Celiac tTG ELISAs were compared with endomysial antibody (EMA) results obtained using a commercially available IFA and a well-characterized clinical population.

### CD population confirmed by EMA

		Positive	Negative	Total
<b>IMMCO</b>	Positive	87	6	93
<b>CELIAC tTG</b>	Negative	1	91	92
<b>IgA ELISA</b>	Total	88	97	185

Positive Percent Agreement: 98.9% (95% CI 92.9% to 99.9%)

Negative Percent Agreement: 93.8% (95% CI 86.5% to 97.5%)

Overall Percent Agreement: 96.2% (95% CI 92.1% to 98.3%)

EMA IgA Positive Celiac Subjects: 88

IgA Deficient Celiac Subjects: 5

Disease Controls: 31

Healthy Normal Subjects: 61

### CD population confirmed by EMA

		Positive	Negative	Total
<b>IMMCO</b>	Positive	72	7	79
<b>CELIAC tTG</b>	Negative	6	191	197
<b>IgG ELISA</b>	Total	78	198	276

Positive Percent Agreement: 92.3% (95% CI 83.4% to 96.8%)

Negative Percent Agreement: 96.5% (95% CI 92.6% to 98.4%)

Overall Percent Agreement: 95.3% (95% CI 91.9% to 97.4%)

EMA Positive Celiac Subjects: 132 (74 EMA IgG Positive)

Disease Controls: 53

Healthy Normal Subjects: 91

C. Clinical Sensitivity and Specificity: Well-characterized CD clinical populations were tested with the Celiac tTG ELISAs producing the following results.

### Celiac Disease

		Positive	Negative	Total
<b>IMMCO</b>	Positive	88	5	93
<b>CELIAC tTG</b>	Negative	5	87	92
<b>IgA ELISA</b>	Total	93	92	185
Sensitivity:		94.6% (95% CI 87.3% to 98.0%)		
Specificity:		94.6% (95% CI 87.2% to 98.0%)		
Percent Agreement:		94.6% (95% CI 90.0% to 97.2%)		
EMA IgA Positive Celiac Subjects: 88				
IgA Deficient Celiac Subjects: 5				
Disease Controls: 31				
Healthy Normal Subjects: 61				

### Non-IgA Deficient CD

		Positive	Negative	Total
<b>IMMCO</b>	Positive	87	0	87
<b>CELIAC tTG</b>	Negative	1	0	1
<b>IgA ELISA</b>	Total	88	0	88
Sensitivity:		98.9% (95% CI 92.9% to 99.9%)		
Specificity:		NA		
Percent Agreement:		98.9% (95% CI 92.9% to 99.9%)		
EMA IgA Positive Celiac Subjects: 88				
Disease Controls: 0				
Healthy Normal Subjects: 0				

### Celiac Disease

		Positive	Negative	Total
<b>IMMCO</b>	Positive	72	7	79
<b>CELIAC tTG</b>	Negative	60	137	197
<b>IgG ELISA</b>	Total	132	144	276
Sensitivity:		54.5% (95% CI 45.7% to 63.2%)		
Specificity:		95.1% (95% CI 89.9% to 97.9%)		
Percent Agreement:		75.7% (95% CI 70.1% to 80.6%)		
CD Subjects Confirmed by EMA: 132 (74 EMA IgG Positive)				
Disease Controls: 53				
Healthy Normal Subjects: 91				

### IgA Deficient CD

		Positive	Negative	Total
<b>IMMCO</b>	Positive	21	0	21
<b>CELIAC tTG</b>	Negative	1	0	1
<b>IgG ELISA</b>	Total	22	0	22
Sensitivity:		95.5% (95% CI 75.1% to 99.8%)		
Specificity:		NA		
Percent Agreement:		95.5% (95% CI 75.1% to 99.8%)		

IgA Deficient CD Subjects: 22

Disease Controls: 0

Healthy Normal Subjects: 0

D. Cross Reactivity: A total of 63 potentially cross-reactive specimens from individuals with other autoimmune disorders or positive for other autoantibodies were tested for tTG antibodies using the Immulisa™ Celiac tTG system.

Condition	n	IgA Positive n (%)	IgG Positive n (%)
Graves Disease	11	0 (0%)	0 (0%)
Hashimoto's Thyroiditis	10	0 (0%)	0 (0%)
ANA positive*	9	0 (0%)	1 (11%)
CCP positive**	10	0 (0%)	0 (0%)
RF positive***	9	0 (0%)	0 (0%)
Ulcerative Colitis	5	0 (0%)	0 (0%)
Crohn's Disease	9	0 (0%)	2 (22%)
Total	63	0 (0%)	3 (5%)

\* Antinuclear Antibodies

\*\* Cyclic Citrullinated Peptides Antibodies

\*\*\* Rheumatoid Factor Antibodies

### Precision

Precision was tested with multiple specimens selected throughout the range of the assay. Three assay runs were performed on different days to determine results between days. An additional run of ten replicates was performed to determine repeatability. Results are summarized below.

Kit	S #	Mean (EU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
			SD (EU/ml)	CV%	SD (EU/ml)	CV%	SD (EU/ml)	CV%
Celiac tTG IgA Assay	1	10.6	1.086	10.2%	1.135	10.7%	1.095	10.2%
	2	20.4	1.224	6.0%	1.476	7.2%	0.956	4.7%
	3	25.0	1.496	6.0%	1.740	6.9%	1.271	5.1%
	4	34.9	2.075	5.9%	2.099	6.1%	2.110	6.0%
	5	48.7	2.634	5.4%	3.661	7.5%	1.090	2.2%
	6	110.6	4.008	3.6%	4.221	3.8%	3.960	3.6%
	7	134.7	7.253	5.4%	8.728	6.5%	5.904	4.4%
	8	183.7	6.161	3.4%	6.606	3.6%	5.615	3.0%
Celiac tTG IgG Assay	1	15.3	1.100	7.2%	1.283	8.3%	0.929	6.1%
	2	24.3	1.727	7.1%	2.101	8.6%	1.355	5.6%
	3	29.8	1.931	6.5%	2.492	8.4%	1.289	4.3%
	4	60.5	3.046	5.0%	3.647	6.0%	2.507	4.1%
	5	99.1	4.720	4.8%	5.687	5.7%	3.756	3.8%
	6	233.7	10.022	4.3%	11.918	5.0%	7.108	3.1%
	7	263.5	13.061	5.0%	14.124	5.4%	11.794	4.4%

### Limit of Detection

The limit of detection (LoD) was determined based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples. LoD for IgA was 4.0 EU/ml. LoD for IgG was 2.9 EU/ml.

### Linearity and Recovery

Linearity and recovery were tested by diluting positive specimens through the assay range in equidistant dilutions and comparing actual vs. expected results. The linear range of the assays were determined be 4.0 (LoD) – 160 EU/ml for IgA and 2.9 (LoD) – 260 EU/ml for IgG. Results are summarized below:

Test Range (EU/ml)	Slope (95% CI)	Y-intercept (95% CI)	R <sup>2</sup>	% recovery (obtained/expected)
IgA				
6.0 to 62.7	0.932 (0.879 to 0.985)	0.926 (-1.240 to 3.091)	0.9943	97.9 to 108.3
3.0 to 157.0	0.947 (0.892 to 1.002)	0.999 (-4.347 to 6.345)	0.9979	94.5 to 107.7
3.3 to 163.6	0.910 (0.824 to 0.957)	3.239 (-1.628 to 8.106)	0.9946	85.0 to 109.5
IgG				
3.0 to 41.9	0.987 (0.970 to 1.004)	0.405 (-0.037 to 0.847)	0.9997	93.5 to 101.2
2.7 to 67.6	0.976 (0.839 to 1.112)	4.303 (-0.854 to 9.459)	0.9801	77.6 to 100.1
2.9 to 260.5	0.800 (0.648 to 0.951)	11.6 (-9.4 to 32.6)	0.9693	85.6 to 122.0

**Interference**

Interference was studied by mixing sera with known tTG antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342  $\mu\text{mol/L}$ ), and Rheumatoid Factor (100 EU/ml).

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1. Guandalini S, Gupta P. Celiac disease- A diagnostic challenge. *Clin Appl Immun Rev* 2002; 2:293-305.
2. Fasano A, et al. Prevalence of celiac disease in at risk and not-at-risk groups in the United States 2003; 163:286-292.
3. Carpenter H, et al. Gastroscopy is incomplete without biopsy: Clinical relevance of distinguishing gastropathy from gastritis. *Gastroenterology* 1995; 108:917-924.
4. Catassi, et al. Risk of non-Hodgkin lymphoma in celiac disease *JAMA* 2002; 287:1413-1419.
5. Kurpa K et al. Diagnosing Mild Enteropathy Celiac Disease: A Randomized, Controlled Clinical Study. *Gastroenterology* 2009; 136: 816-823.
6. Hill ID, et al. Celiac disease: Working group report of the first world congress of pediatric gastroenterology, hepatology, and nutrition. *J. Pediatr Gastroenterol Nutr* 2002; 35:S78-S88.
7. Rosario D, et al. Further studies of anti-endomysium and anti-gliadin antibodies in patients with suspected celiac disease. *J. Pediatr Gastroenterol Nutr* 1998; 27:191-195.
8. Sárdy M, et al. Tissue transglutaminase ELISA positivity in autoimmune disease independent of gluten-sensitive disease. *Clin Chim Acta.* 2007;376:126-35;
9. Song KS et al. Tissue transglutaminase autoantibodies in patients with IgM rheumatoid factors. *Yonsei Med J.* 2004;45:960-2.
10. Bizzaro N et al. Low specificity of anti-tissue transglutaminase antibodies in patients with primary biliary cirrhosis. *J Clin Lab Anal.* 2006;20:184-9.
11. Villalta D et al. False positive reactions for IgA and IgG anti-tissue transglutaminase antibodies in liver cirrhosis are common and method-dependent. *Clin Chim Acta.* 2005;356:102-09.
12. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
13. Korponay-Szabo IR, Dahlbom I, Laurila K, et al. Elevation of IgG antibodies against tissue transglutaminase as a diagnostic tool for coeliac disease in selective IgA deficiency. *Gut.* 2003;52:1567-1571.
14. Sulkanen S, Haltunen T, Laurila K, et al. Tissue Transglutaminase autoantibody enzyme-linked immunosorbent assay in detecting celiac disease. *Gastroenterology*; 1998, 115:1322-1328.
15. Dietrich W, Laag E, Schopper H, et al. Autoantibodies to tissue transglutaminase as predictors of celiac disease; *Gastroenterology* 1998, 115:1317-1321.


*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**  
**60 Pineview Drive**  
**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**  
**Fax: (716) 691-0466**  
**Toll Free USA/Canada: 1-800-537-TEST**  
**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands  
Tel (+31) 345 8570, Fax (+31) 346 7299  
[www.emergogroup.com](http://www.emergogroup.com)

DEC2011

Document No. PI5144