For technical assistance please contact:

**IMMCO Diagnostics, Inc.**  
60 Pineview Drive  
Buffalo, NY 14228-2120  
Telephone: (716) 691-0091  
Fax: (716) 691-0466  
Toll Free USA/Canada: 1-800-537-TEST  
E-Mail: info@immcodiagnostics.com

or your local product distributor

EU Authorized Representative/ Εξουσιοδοτημένος αντιπρόσωπος της / Representante Autorizado/ Autorisierter Repräsentant/ Représentant Autorisé /Rappresentante Autorizzato / Representante Autorizado  
EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands  
Tel (+31) 345 8570, Fax (+31) 346 7299  
www.emergogroup.com

MAR2010

Document No. PI5159

# ImmunoLISA™ Celiac G+

## Gliadin Peptide Antibody ELISA

**IVD** For *in vitro* diagnostic use

### PRODUCT INSERT

<b>REF</b>	5159A	Celiac G+ IgA ELISA	96 Determinations
<b>REF</b>	5159G	Celiac G+ IgG ELISA	96 Determinations

### INTENDED USE

Enzyme linked immunosorbent assay (ELISA) for the qualitative and semi-quantitative detection of IgA or IgG antibodies to gliadin in human serum to aid in the diagnosis of celiac disease in conjunction with other laboratory and clinical findings.

### SUMMARY AND EXPLANATION

Celiac Disease (CD) is an autoimmune gastrointestinal disorder that may occur in genetically susceptible individuals triggered by the ingestion of gluten-containing grains such as wheat, barley and rye. CD is characterized by malabsorption resulting from inflammatory injury to the small intestinal mucosa and, when prolonged, can cause malnutrition. The classical symptoms of CD include diarrhea, weight loss and malnutrition. Only a small percentage of patients with CD, however, presents with classical symptoms.<sup>1</sup> Consequently, the clinical spectrum of CD has grown much broader than in the past to include patients that do not present with classical symptoms. It is not uncommon for the initial symptoms to be non-gastrointestinal or for gastrointestinal symptoms, if present, to be mild or intermittent. The need to examine a wider range of clinical presentation has led to greater numbers of individuals diagnosed with CD later in life than ever before. The advent of serological methods for the detection of antibodies to gliadin, endomysium and tissue transglutaminase have enabled large-scale screening studies for CD both in Europe and the United States. These studies suggest that CD is far more prevalent than thought previously. Recent serological studies demonstrate incidence of CD between one in 100 and one in 500.<sup>2,3</sup> Prevalence of CD is much higher in first and second-degree relatives of patients with CD. CD has been associated with many other autoimmune disorders such as type 1 diabetes, thyroid autoimmunity and other autoimmune disorders.

The most common serologic tests for the screening of CD are the immunofluorescence method of detecting endomysial antibodies (EMA) and ELISA methods of detecting antibodies to tissue transglutaminase (tTG) and gliadin.<sup>4-8</sup> Because of the limited sensitivity and specificity of gliadin antibody assays this method has largely been replaced with EMA and tTG.

The limitation of the EMA and tTG immunoassays is that both immunoassays detect IgA antibodies, hindering identification of IgA deficient patients with CD.<sup>9-11</sup> IgG gliadin antibody tests are important towards the diagnosis of CD in patients who are IgA deficient.<sup>10</sup> Studies show that 1-2% of the general population is IgA

deficient and that the incidence of CD in IgA deficient subjects is significant, hence the need for specific tests. In addition, some studies question the specificity of tTG assays and the sensitivity of EMA.<sup>12-16</sup> Gliadin peptide based assays have shown promising results in addressing the limitations of current assays.<sup>17-23</sup> This next generation assay for the detection of gliadin antibodies, the Celiac G+ ELISA, incorporates a proprietary gliadin peptide technology that provides sensitive and specific detection of IgA or IgG isotype antibodies with high degree of reliability.

## PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with deamidated gliadin peptides followed by a blocking step to reduce non-specific protein binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the gliadin antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgA or IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml) and reported as positive or negative.

## REAGENTS

### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

### Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.<sup>24</sup>

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other

**REFERENCES/ΠΑΡΑΠΟΜΠΕΣ/REFERENCIAS/REFERENZEN/  
RÉFÉRENCES/RIFERIMENTI/REFERÊNCIAS**

1. Guandalini S et al. *Clin Appl Immun Rev*. 2002;2:293-305.
2. Fasano A, et al. *Arch Intern Med*. 2003; 163:286-292.
3. Lohi S et al. *Aliment Pharmacol Ther*. 2007; 26:1217-25.
4. Chorzelski TP et al. *Br J Dermatol*. 1984;111:395-402.
5. Pacht A et al. *Isr J Med Sci*. 1995; 31:218-20.
6. Grodzinsky E et al. *Acta Paediatr*. 2008; 97:972-6.
7. Tucker NT et al. *J Pediatr*. 1988; 113:286-289.
8. Schuppan D et al. *Eur J Gastroenterol Hepatol*. 2001; 13: 635-7.
9. Korponay-Szabo IR et al. *Gut*. 2003; 52:1567-1571.
10. Kumar V et al. *Clinical Diagnostic Lab Immunology*. 2002; 9:1295-1300.
11. McGowan KE et al. *Clin Chem*. 2008; 54:1203-9.
12. Song KS, Choi JR. *Yonsei Med J*. 2004; 45:960-2.
13. Luft LM et al. *J Rheumatol*. 2003; 30:2613-9.
14. Sárdy M et al. *Clin Chim Acta*. 2007; 376:126-35.
15. Bizzaro N et al. *J Clin Lab Anal*. 2006; 20:184-9.
16. Villalta D et al. *Clin Chim Acta*. 2005; 356:102-9.
17. Mothes T. *Adv Clin Chem*. 2007; 44:35-63.
18. Kaukinen K et al. *Scand J Gastroenterol*. 2007; 42:1428-33.
19. Schwertz E et al. *Clin Chem*. 2004; 50: 2370-5.
20. Rashtak S et al. *Clin Gastroenterol Hepatol*. 2008; 6:426-32
21. Volta U et al. *Dig Dis Sci*. 2007; 3:1582-8.
22. Ankelo M et al. *Clin Exp Immunol*. 2007; 150: 285-293.
23. Bansal A et al. *Annals NY Acad Sci*. 2009.
24. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).

sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

**Materials provided**

ImmuliS<sup>™</sup> Celiac G+ IgA ELISA REF 5159A  
 ImmuliS<sup>™</sup> Celiac G+ IgG ELISA REF 5159G

Kits contains sufficient reagents to perform 96 determinations.

12 x 8	MICROPLATE CG+	Microplate with individual breakaway microwells. Coated with deamidated gliadin peptides. Ready for use.
1 x 1.75 ml	CONTROL+ CG+A	Ready to use <b>Positive Control</b> (red cap) for <span style="border: 1px solid black; padding: 2px;">REF</span> 5159A. Contains human serum positive for gliadin IgA antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL+ CG+G	Ready to use <b>Positive Control</b> (red cap) for <span style="border: 1px solid black; padding: 2px;">REF</span> 5159G. Contains human serum positive for gliadin IgG antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL-	Ready to use <b>Negative Control</b> (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A CG+A	Ready to use set of 5 <b>Calibrators</b> for <span style="border: 1px solid black; padding: 2px;">REF</span> 5159A. Calibrator A (green cap) 320 EU/ml, Calibrator B (violet cap) 160 EU/ml, Calibrator C (blue cap) 80 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing gliadin IgA antibodies. Concentrations in EU/ml are printed on the labels.
	CALIBRATOR B CG+A	
	CALIBRATOR C CG+A	
	CALIBRATOR D CG+A	
	CALIBRATOR E CG+A	
5 x 1.75 ml	CALIBRATOR A CG+G	Ready to use set of 5 <b>Calibrators</b> for <span style="border: 1px solid black; padding: 2px;">REF</span> 5159G. Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing gliadin IgG antibodies. Concentrations in EU/ml are printed on the labels.
	CALIBRATOR B CG+G	
	CALIBRATOR C CG+G	
	CALIBRATOR D CG+G	
	CALIBRATOR E CG+G	
1 x 15 ml	IgA-CONJ HRP	HRP goat anti-human <b>IgA Conjugate</b> for <span style="border: 1px solid black; padding: 2px;">REF</span> 5159A. Ready for use. Color coded pink.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human <b>IgG Conjugate</b> for <span style="border: 1px solid black; padding: 2px;">REF</span> 5159G. Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. <b>Protect from light.</b>
1 x 15 ml	STOP H2SO4	Stop Solution. Ready for use.
2 x vials	BUF WASH	Powder Wash Buffer. <b>Reconstitute to one liter each.</b>
1 x		Protocol Sheets

**Optional Components**

1 x 60ml BUF|WASH Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

## Symbols used on labels

	Lot number
	Catalog number
	In vitro diagnostic use
	Use by
	Storage temperature
	Read instructions before use
	Number of tests
	Manufacturer

## Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

## SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2° - 8°C for no longer than one week. For longer storage, serum specimens should be frozen. It is recommended that frozen specimens be tested within one year. Avoid repeated freezing and thawing of samples.

## PROCEDURE

### Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.

Série do teste (EU/ml)	Inclinação (95% CI)	Intercepção Y (95% CI)	R <sup>2</sup>	% de recuperação (obtida/prevista)
<b>IgA</b>				
7,2 - 70,3	0,971 (0,907 - 1,035)	0,468 (-2,415 - 3,349)	0,9958	97,3 - 106,9
4,4 - 136	0,987 (0,942 - 1,033)	2,203 (-1,518 - 5,924)	0,9979	88,9 - 103,3
1,1 - 309,6	0,993 (0,909 - 1,076)	-4,609 (-18,092 - 11,808)	0,9946	94,4 - 114,8
<b>IgG</b>				
5,2 - 75,8	1,011 (0,975 - 1,047)	-0,499 (-2,170 - 1,171)	0,9987	97,3 - 105,7
4,8 - 93,9	1,006 (0,946 - 1,067)	1,161 (-2,217 - 4,538)	0,9963	90,9 - 102
1,2 - 169,8	1,034 (0,992 - 1,076)	0,9027 (-3,196 - 5,002)	0,9971	86,0 - 102,3

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- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

### Test Method

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.
- or**
- For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

	Qualitative			Semi-Quantitative			
A	Blank	S5	Etc.	A	Blank	S1	Etc.
B	-	S6		B	-	S2	
C	+	S7		C	+	S3	
D	Cal D	S8		D	Cal A	S4	
E	S1	S9		E	Cal B	S5	
F	S2	S10		F	Cal C	S6	
G	S3	S11		G	Cal D	S7	
H	S4	S12		H	Cal E	S8	
	1	2	3		1	2	3

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500µl** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** ( $\pm$  5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To

blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.

- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <5 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

## RESULTS

### Calculations

The concentrations of the patient samples can be determined by either of two methods:

#### 1. QUALITATIVE DETERMINATION

Abs. of Test Sample

----- X EU/ml of Calibrator D = EU/ml Test Sample

Abs. of Calibrator D

It is recommended that qualitative results be reported as "positive" or "negative." Sample results greater than or equal to Calibrator D are considered positive.

#### 2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve

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according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with EU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods, such as assays for detection of EMA and/or tTG antibodies.

#### Interpretation

The following serves only as a guide in the interpretation of the laboratory results. These values were determined by testing 64 normal human sera and 50 non-celiac disease controls. The values depicted below are the mean of the normal subjects plus 2SD. Each laboratory must determine its own normal values.

anti-Gliadin Ab value	Interpretation
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

#### Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

#### LIMITATIONS OF PROCEDURE

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

#### EXPECTED VALUES

Test results in a normal population are usually expected to be negative. However, as the incidence of CD in normal population is about 1%, some apparently healthy, asymptomatic individuals may test positive on the Celiac G+ ELISA.

#### PERFORMANCE CHARACTERISTICS

The utility of the ImmuLisa™ Celiac G+ ELISAs for the detection of gliadin antibodies was evaluated by testing well-characterized EMA positive serum specimens from suspected CD subjects alongside disease controls and “normal” human sera. These specimens were also tested on commercially available test kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

A. Celiac G+ vs. other Gliadin Peptide Antibody kit:

**Other Gliadin IgA ELISA**

		Positive	Negative	Total
<b>IMMCO</b>	Positive	54	27	81
<b>CELIAC G+</b>	Negative	7	68	75
<b>IgA ELISA</b>	Total	61	95	156

Positive Percent Agreement: 88.5% (95% CI 77.2% - 94.9%)  
 Negative Percent Agreement: 71.6% (95% CI 61.3% - 80.1%)  
 Overall Percent Agreement: 78.2% (95% CI 70.8% - 84.3%)

**Other Gliadin IgG ELISA**

		Positive	Negative	Total
<b>IMMCO</b>	Positive	88	9	97
<b>CELIAC G+</b>	Negative	7	109	116
<b>IgG ELISA</b>	Total	95	118	213

Positive Percent Agreement: 92.6% (95% CI 84.9% - 96.7%)  
 Negative Percent Agreement: 92.4% (95% CI 85.6% - 96.2%)  
 Overall Percent Agreement: 92.5% (95% CI 87.9% - 95.5%)

B. Celiac G+ vs. EMA: Results obtained with the Celiac G+ ELISAs were compared with endomysial antibody (EMA) EMA titer obtained using a commercially available IFA kit and a well-characterized clinical population. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

**CD population confirmed by EMA**

		Positive	Negative	Total
<b>IMMCO</b>	Positive	78	5	83
<b>CELIAC G+</b>	Negative	18	55	73
<b>IgA ELISA</b>	Total	96	60	156

Estimated Sensitivity: 81.3% (95% CI 71.7% - 88.2%)  
 Estimated Specificity: 91.7% (95% CI 80.9% - 96.9%)  
 Agreement: 85.3% (95% CI 78.5% - 90.2%)

EMA Positive Celiac Subjects: 96  
 Disease Controls: 7  
 Healthy Normal Subjects: 53

**CD population confirmed by EMA**

		Positive	Negative	Total
<b>IMMCO</b>	Positive	94	3	97
<b>CELIAC G+</b>	Negative	11	105	116
<b>IgG ELISA</b>	Total	105	108	213

Estimated Sensitivity: 89.5% (95% CI 81.6% - 94.4%)  
 Estimated Specificity: 97.2% (95% CI 91.5% - 99.3%)  
 Agreement: 93.4% (95% CI 89.0% - 96.2%)

EMA Positive Celiac Subjects: 105  
 Disease Controls: 20  
 Healthy Normal Subjects: 88

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C. Cross Reactivity: A total of 50 potentially cross-reactive specimens from individuals with other autoimmune disorders or positive for other autoantibodies were tested for Gliadin antibodies using the Immulisa™ Celiac G+ system.

Condition	n	IgA Positive n (%)	IgG Positive n (%)
Graves' Disease	11	0 (0%)	0 (0%)
Hashimoto's Thyroiditis	11	1 (9%)	0 (0%)
ANA positive*	9	0 (0%)	0 (0%)
CCP positive**	10	0 (0%)	2 (20%)
RF positive***	9	1 (11%)	0 (0%)
Total	50	2 (4%)	2 (4%)

\* Antinuclear Antibodies

\*\* Cyclic Citrullinated Peptides Antibodies

\*\*\* Rheumatoid Factor Antibodies

### Precision

Precision was tested with multiple specimens selected throughout the range of the assay. Three assay runs were performed on different days to determine results between days. An additional run of six replicates was performed to determine repeatability. Results are summarized below.

Kit	S #	Mean (EU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (EU/ml)	CV%	SD (EU/ml)	CV%	SD (EU/ml)	CV%	
Celiac G+ IgA Assay	1	15.38	1.713	11.1%	2.474	16.2%	0.572	3.7%
	2	22.09	1.649	7.5%	2.011	9.1%	1.392	6.3%
	3	27.23	1.792	6.6%	1.932	7.1%	1.807	6.7%
	4	37.08	2.565	6.9%	3.043	8.2%	2.277	6.2%
	5	79.94	6.027	7.5%	8.438	10.7%	2.686	3.3%
	6	180.68	5.734	3.2%	7.852	4.3%	3.205	1.8%
	7	186.77	7.926	4.2%	10.837	5.9%	3.872	2.1%
	8	341.97	9.786	2.9%	10.389	3.0%	10.137	3.0%
Celiac G+ IgG Assay	1	13.97	0.827	5.9%	0.884	6.4%	0.757	5.3%
	2	23.10	1.274	5.5%	0.327	1.4%	1.821	7.8%
	3	33.81	2.071	6.1%	2.552	7.6%	1.696	5.0%
	4	91.83	3.472	3.8%	4.504	5.0%	1.723	1.9%
	5	116.49	4.360	3.7%	5.813	5.1%	1.349	1.1%
	6	120.44	3.008	2.5%	2.713	2.3%	2.770	2.3%
	7	179.74	4.661	2.6%	5.461	3.1%	3.998	2.2%

## Limit of Detection

The limit of detection (LoD) was determined based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples. LoD for IgA was 8.9 EU/ml. LoD for IgG was 2.2 EU/ml.

## Linearity and Recovery

Linearity and recovery were tested by diluting positive specimens through the assay range in equidistant dilutions and comparing actual vs. expected results. The linear range of the assays were determined be 8.9 (LoD) – 320 EU/ml for IgA and 2.2 (LoD) – 160 EU/ml for IgG. Results are summarized below:

Test Range (EU/ml)	Slope (95% CI)	Y-intercept (95% CI)	R <sup>2</sup>	% recovery (obtained/expected)
<b>IgA</b>				
7.2 to 70.3	0.971 (0.907 to 1.035)	0.468 (-2.415 to 3.349)	0.9958	97.3 to 106.9
4.4 to 136	0.987 (0.942 to 1.033)	2.203 (-1.518 to 5.924)	0.9979	88.9 to 103.3
1.1 to 309.6	0.993 (0.909 to 1.076)	-4.609 (-18.092 to 11.808)	0.9946	94.4 to 114.8
<b>IgG</b>				
5.2 to 75.8	1.011 (0.975 to 1.047)	-0.499 (-2.170 to 1.171)	0.9987	97.3 to 105.7
4.8 to 93.9	1.006 (0.946 to 1.067)	1.161 (-2.217 to 4.538)	0.9963	90.9 to 102
1.2 to 169.8	1.034 (0.992 to 1.076)	0.9027 (-3.196 to 5.002)	0.9971	86.0 to 102.3

## Interference

Interference was studied by mixing sera with known gliadin antibody levels with potentially interfering serum samples and studying deviation from expected results. Results for IgA and IgG isotypes appear below.

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