

SUMMARY OF PROCEDURE

1. Prepare a 1:101 dilution of patient samples in Sample Diluent. Mix well. Note that Calibrator and Controls are ready-to-use.
2. Add 100 µl of Calibrator, Controls and diluted samples into the antigen wells. Reserve one well for reagent blank (100 µl of Sample Diluent).
3. Incubate at room temperature (18-30°C) for 30 ± 5 min.
4. Discard contents of the wells. Wash the wells 3 times with Wash Solution.
5. Add 100 µl of Conjugate to each well.
6. Incubate at room temperature (18-30°C) for 30 ± 5 min.
7. Wash the wells as in #4 above.
8. Add 100 µl Substrate Solution to each well.
9. Incubate at room temperature (18-30°C) for 30 ± 5 min.
10. Add 100 µl Stop Solution to each well.
11. Read the absorbances at 450/600-630 nm.

INTENDED USE

For the semi-quantitative detection of IgG antibodies to measles virus in human serum by indirect enzyme immunoassay to aid in the assessment of the patient's immunological response to measles virus. Additionally, it is intended to aid in the qualitative determination of the immune status of individuals and, when evaluating paired sera, as an aid in the diagnosis of measles infection.

SUMMARY AND EXPLANATION

Measles (Rubeola) is a highly contagious infection caused by an RNA myxovirus. The incubation period is 10 to 11 days and the infection is characterized by fever, myalgias, nonproductive cough, conjunctivitis and exanthem and enanthem (Koplik's spots). The rash of rubeola almost always begins on the face and then spreads to the trunk and extremities. Typically, the illness crests on day three of fever and the temperature falls to normal on day seven (1,2).

Prior to the introduction of vaccines, measles was an inevitable disease of childhood. However, since the introduction of the measles vaccine in 1963, the incidence of this disease has dropped dramatically and physicians have become less familiar with the disease (3). Diagnosis of the disease can become further complicated by the emergence of atypical forms of measles. These atypical forms generally occur in recipients of inactivated measles vaccine who were immunized with this vaccine in 1963-1967 and were subsequently exposed to the natural disease. The atypical form of measles may be severe and is often confused with Rocky Mountain spotted fever (1). In addition, acute measles infection can be complicated by secondary infections of the lower respiratory tract and ear. Additional complications such as encephalomyelitis occur in about 0.1% of patients (3). Measles infection during pregnancy has been associated with an increased risk of miscarriages or premature delivery (4). Persistent measles infection has also been suspected in chronic autoimmune disease. Patients with multiple sclerosis, as well as those with systemic lupus erythematosus, have consistently been found to have elevated levels of antibodies to measles virus. This does not necessarily indicate a role of the virus in the etiology of these diseases, but may reflect polyclonal activation of B cells known to occur during the courses of these diseases (3).

The presence of specific antibodies in a single serum specimen indicates past measles infection or vaccination. Demonstration of a significant increase in antibody titers in a serum pair taken at a 7-14 day interval is the basis for diagnosis of acute infection (3).

The traditional methods of antibody detection such as hemagglutination inhibition (HI), and neutralization (Nt) have been replaced by the enzyme immunoassay (EIA) which is more sensitive, equally specific and less labor intensive (5,6).

The Diamedix Immun simplicity® Measles IgG Test Kit is an EIA procedure intended for the semi-quantitative detection of antibodies to measles virus antigen. The results are objective and reported in ELISA units (EU/ml), standardized against in-house reference materials.

PRINCIPLE OF THE PROCEDURE

Diluted samples are incubated with measles virus antigen bound to the solid surface of a microtiter well. If IgG antibodies against measles are present in the samples, they will bind to the antigen forming antigen-antibody complexes. Residual sample is eliminated by aspirating and washing. Conjugate (horseradish peroxidase-labeled anti-human IgG) is added and will bind to these complexes. Unbound conjugate is removed by aspiration and washing. Substrate is then added and incubated. In the presence of bound enzyme the substrate is converted to an end product. The absorbance of this end product can be read spectrophotometrically at 450 nm (reference 600-630 nm) and is directly proportional to the concentration of IgG antibodies to measles antigen present in the sample.

REAGENTS

Each Is-Measles IgG Test Kit contains reagents for 96 tests.

Antigen Wells	Twelve, 8-well microwell breakapart strips, color-coded dark blue, coated with Measles virus antigen (partially purified extract of Vero cells infected with the Edmonston strain of measles virus).
Calibrator	Two vials with blue caps containing 1.8 ml of pre-diluted human serum, highly reactive for measles IgG antibodies, 0.2% sodium azide and Proclin® 300, 90 ppm active ingredient. Assigned EU/ml value printed on label.
Positive Control	One vial with white cap containing 1.8 ml of pre-diluted human serum, moderately reactive for measles IgG antibodies, 0.2% sodium azide and Proclin® 300, 90 ppm active ingredient. Assigned EU/ml range printed on label. Note that controls and Calibrators are produced from different serum lots.
Negative Control	One vial with black cap containing 1.8 ml of pre-diluted human serum, non-reactive for measles IgG antibodies, 0.2% sodium azide and Proclin® 300, 90 ppm active ingredient.
Sample A Diluent	One bottle with blue cap containing 60 ml Phosphate buffer with protein stabilizers. Contains 0.2% sodium azide and Proclin® 300, 90 ppm active ingredient. Color-coded blue.
Wash S Concentrate (20X)	Two bottles with clear caps containing 50 ml of Phosphate buffered saline with Proclin® 300, 15 ppm active ingredient. Color-coded light blue/green. Each bottle is sufficient to make 1050 ml liter of wash solution.
Conjugate	One bottle with red cap containing 25 ml goat anti-human immunoglobulin G labeled with horseradish peroxidase. Also includes protein stabilizers and preservatives. Color-coded pink.
Substrate HRP	One amber bottle with brown cap containing 25 ml buffered TMB solution (3, 3', 5, 5' tetramethylbenzidine dimethylsulfoxide).
Stop M Solution	One bottle with white cap containing 30 ml of 1 N Phosphoric and 1N Hydrochloric acids. CAUTION: Acids are corrosive. Avoid contact with skin or eyes. If contact is made, flush area with copious amounts of water. See Precautions section.

Store these reagents at 2 to 8°C.

OTHER MATERIALS REQUIRED

Manual Users:

- Wash bottle or automated microplate washer.
- Pipettors capable of dispensing appropriate volumes.
- Timer.
- One liter graduated cylinder.
- One liter wash solution reservoir.
- Deionized or distilled water.
- Absorbent toweling.
- Tubes or microwell plate for serum dilution.
- Reader capable of reading absorbance at 450 nm, reference at 600-630 nm.

Automated EIA Processor Users:

- One liter graduated cylinder.
- Deionized or distilled water.
- Pre-dilution cups, strips or plates.
- ProbeClean™ Concentrate, or tip washing detergent solution, if applicable.

PRECAUTIONS

REAGENTS: For *in vitro* Diagnostic Use.

1. Handle samples, Calibrator, controls and the materials that contact them as potential biohazards. Each donor unit in the standards and controls has been found negative for Hepatitis B surface antigen and HIV-1 antibodies by FDA-approved third generation tests. However, because no method can offer complete assurance that HIV-1, Hepatitis B virus, or other infectious agents are absent, these materials should be handled at the Biosafety Level 2 as recommended for any potentially infectious serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories", 1993.
2. Never pipette by mouth.
3. Avoid contact with open skin and mucous membranes.
4. Certain of the test reagents contain Proclin[®] 300 as a preservative. When disposing of reagents containing Proclin[®] 300, flush drains with copious amounts of water to dilute the components below active levels.
5. Reagents containing Sodium Azide:
 - (a) **CAUTION:** Some reagents in this kit contain Sodium Azide as preservative. Sodium Azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts", in the Manual Guide – Safety Management No. CDC-22, issued by the Centers for Disease Control and Prevention, Atlanta, GA, 1976.

European Communities Hazardous Substance Risk Phrases (Regulation (EC) No 1272/2008)

H300 –Fatal if swallowed.
H310 – Fatal if contact with skin.
EUH032 – Contact with acids liberates very toxic gas.
H410 – Very toxic to aquatic life with long lasting effect.
P264 – Wash all exposed external body areas thoroughly after handling.
P302+P352 – IF ON SKIN: Wash with plenty of water and soap.
P301+P310/P330 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Rinse mouth.
P270 – Do not eat, drink or smoke when using this product.
P501 – Dispose of contents/container as hazardous waste.
P391 – Collect spillage.
P273 – Avoid release to the environment. Refer to special instructions/ Safety Data Sheet.
 - (b) Sodium Azide inhibits horseradish peroxidase activity. Care must be taken to ensure that azide is not carried over from other reagents into conjugate and substrate steps.
6. Avoid contamination of the TMB substrate solution with conjugate or other oxidants which will cause the solution to change color prematurely.

ADDITIONAL PRECAUTIONS:

1. Do not interchange reagents from different reagent lots except for Sample Diluent, Wash Concentrate, Substrate and Stop Solution.
2. Do not use reagents beyond their expiration date. Expiration dates are printed on the reagent labels.
3. Store unused reagents at 2 to 8°C.
4. Incubations above or below the recommended temperatures or times may give erroneous results.
5. The EIA method is a very sensitive technique. Maintain consistent pipetting technique, incubation times, and temperature conditions throughout the test procedure. Cross contamination between reagents can invalidate the test.
6. Antigen coated microwells should be stored with the desiccant in the resealable bag provided and returned to the refrigerator immediately after use.
7. (*Manual Procedure Only*) The washing procedure is very important and requires special attention. (Please refer to the Procedure section).

NOTE: *Improperly washed wells may give erroneous results.*
8. The concentration of anti-Rubeola (measles) IgG in a given specimen determined from assays from different manufacturers can vary due to differences in assay methods and reagents.

SPECIMEN COLLECTION

Whole blood should be collected by accepted medical techniques. Separated serum should remain at 22°C for no longer than 8 hours. If assays are not

completed within 8 hours, serum should be refrigerated (2 to 8°C). If assays are not completed within 48 hours, or the separated sample is to be stored beyond 48 hours, samples are to be frozen at –20°C. Prior to testing, bring frozen sera to room temperature slowly and mix gently, avoiding foam formation. Specimens containing visible particulate matter should be clarified by centrifugation before testing. Grossly contaminated, hemolyzed, lipemic, or icteric specimens should not be used.

For the diagnosis of acute measles infection, the acute-phase specimen should be drawn as soon after onset as possible, preferably within the first 7 days. The convalescent-phase specimen should be drawn 10 or more days after the acute-phase specimen.

CAUTION: Serum samples must not be heat-inactivated prior to use.

PROCEDURE

Allow all test components and patient samples to warm to room temperature before use. Invert reagent bottles gently several times before use. Return promptly to the refrigerator after use.

Prepare Wash Solution by adding 50 ml of Wash Concentrate (20X) to one liter of deionized or distilled H₂O.

MANUAL USERS:

The Calibrator and Controls are provided ready to use: DO NOT DILUTE FURTHER.

1. Prepare 1:101 dilutions of the patient samples in Sample Diluent. (e.g., by addition of 2 μ l sample to 200 μ l Sample Diluent or 5 μ l sample to 500 μ l Sample Diluent).
2. Mix sample dilutions gently by withdrawing and expelling in a pipette tip 2 or 3 times or by vortex mixing for 2 or 3 seconds. Transfer 100 μ l of Calibrator, controls and diluted patient samples, to the antigen wells. Avoid formation of bubbles when transferring diluted samples.

NOTE: *Include one well which contains 100 μ l of Sample Diluent only as the reagent blank. This will ultimately be used to "zero" the photometer before reading the test results.*
3. Allow the wells to incubate at room temperature (18 to 30°C) for 30 \pm 5 minutes.
4. Aspirate or discard the contents of the wells. Remove any excess moisture in the wells by tapping on paper toweling. Wash the wells by rinsing 3 times with Wash Solution. Remove excess moisture from the wells after washing. When using an automated washer, follow the manufacturer's instructions.
5. Place 100 μ l of Conjugate into each well, avoiding bubble formation.
6. Allow the wells to incubate uncovered at room temperature (18 to 30°C) for 30 \pm 5 minutes.
7. Wash the wells as described in Step 4 above.
8. Place 100 μ l of Substrate into each well, avoiding bubble formation.
9. Allow the wells to incubate uncovered at room temperature (18 to 30°C) for 30 \pm 5 minutes.
10. Place 100 μ l of Stop Solution into each well, avoiding bubble formation.
11. Read the absorbance of the wells at 450 nm using a reference wavelength of 600-630 nm. The plate should be read within 60 minutes of adding Stop Solution.

AUTOMATED EIA PROCESSOR USERS:

If using an Automated EIA Processor, refer to the Operator's Manual for the test setup and procedures.

NOTE: *Automated EIA Processor users must validate their equipment to demonstrate that the results obtained are equivalent to those obtained using manual assay.*

QUALITY CONTROL

1. The Positive and Negative Controls must be included in each test run.
2. The absorbance of the Blank must be < 0.2.
3. The Positive Control must be within its assigned range.
4. The Negative Control must be < 15.0 EU/ml.

If any of these criteria are not met, the results are invalid and the test should be repeated.

Notes: The Negative and Positive Controls are intended to monitor substantial reagent failure. The controls will not control all parts of the procedure such as technical dilution of patient specimens. The Positive Control will not ensure precision at the assay cut-off. Users may wish to

establish an in-house control having a quantitative value determined by replicate testing, at or near the cut-off to monitor the precision of the assay cut-off.

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations.

RESULTS

Single Point Calibration

The Diamedix *Is*-Measles IgG Test Kit has been developed using a single point calibrator. The use of a single point calibrator is possible since the test system has been shown to be linear from the cut-off up to the Calibrator value (100 EU/ml). This linear response has been validated both manually and using MAGO[®] by testing serial dilutions of the Calibrator. R² values obtained for both manual and MAGO Plus were greater than 0.99. If users wish to check the linear response they may do so by testing serial two-fold dilutions of the Calibrator.

1. Calculation

Determine the EU/ml (ELISA Units/ml) for each patient specimen or control using the following formula:

$$\frac{\text{EU/ml of Calibrator}}{\text{Absorbance of Calibrator}} \times \text{Absorbance of sample} = \text{EU/ml of sample}$$

Patient values which produce absorbance values above the Calibrator absorbance may be reported as "greater than 100 EU/ml". Alternatively such samples may be pre-diluted using Sample Diluent and re-assayed. Several dilutions (for example 1/5, 1/10 and 1/20) of the pre-diluted sample may be re-assayed simultaneously. Select the dilution that has an absorbance reading about 50% of the absorbance reading of the Calibrator; calculate the EU/ml for this dilution, and multiply by the dilution factor to obtain estimated values.

Index values can be calculated by dividing the EU/ml values obtained by 20 (the positive cut-off value).

An Automated EIA Processor (e.g. MAGO[®] Plus Automated EIA Processor) will calculate and print results automatically using the above formula.

2. Interpretation

EU/ml	Index Value	Interpretation
< 15.0 EU/ml	< 0.75	Nonreactive (Negative) for anti-Rubeola IgG: presumed non-immune to measles virus
15.0 – 19.9 EU/ml	0.75 – 0.99	Equivocal*
≥ 20.0 EU/ml	≥ 1.0	Reactive (Positive) for anti-Rubeola IgG: presumed immune to measles virus

* When equivocal results are obtained, another specimen should be collected ten to fourteen days later and tested in parallel with the initial specimen. If the second sample is also equivocal, the patient is negative for primary or recent infection, and equivocal for antibody status. If the second sample shows a significant increase in antibody level, the patient may be considered to have a primary infection. The conversion of an individual patient's serum from negative to positive for antibodies to the infectious agent in question, is defined as seroconversion, and indicates active or recent infection.

3. Reporting Results

When the EU/ml value is reported for a single specimen the following statement should be included: "The following results were obtained with the Diamedix Immunosimplicity[®] Measles IgG EIA test system. The magnitude of the measured result, above the cut-off, is not indicative of the total amount of antibody present. The magnitude of the reported IgG level cannot be correlated to an endpoint titer".

When the assay is used semi-quantitatively, the following statement should be included when reporting results: "Timing of specimen collection for paired sera may be critical. In some patients, antibody titers may rise to significant levels and fall again to lower or undetectable levels within a month. Other patients may not develop significant antibody levels. Culture results, serology and antigen detection methods should all be appropriately used along with clinical findings for diagnosis".

4. Paired Sera

To determine a significant difference between acute/convalescent serum pairs, both specimens should be run within the same assay. In addition, paired sera should be evaluated within the linear range of the assay. The upper limit of the linear range has been set at 100 EU/ml. In-house studies performed manually and using the MAGO[®] have shown that a 2.1-fold to a 4.4-fold (mean 3.2-fold ± 2 SD) increase in Index Ratio (convalescent serum Index value / acute serum Index value) corresponds to a four-fold increase in measles IgG antibody level. An Index Ratio in the range of 1.5 to 2.1 indicates an equivocal status for the paired sample Index Ratio. In this case, paired samples can be retested or additional samples collected if necessary. If paired sera controls are desired, it is recommended that a four-fold dilution of the Calibrator, or other known positive sample of Calibrator strength, is made first in sample diluent and then diluted according to assay procedures. The undilute and 4-fold diluted material will provide a simulated serum pair. The four-fold dilution Index ratio is compared against the established range (see above).

CUT-OFF ESTABLISHMENT

The Diamedix *Is*-Measles IgG Test Kit cut-off was established to optimally differentiate those individuals with, from those without, immunological experience to measles. The optimal cut-off value was determined by statistical analysis (mean plus 3 SD) of serum samples non-reactive (negative) for measles antibodies. This cut-off value was further verified using the principles from Receiver-Operating Characteristic (ROC) Curves on the two hundred and four samples from normal blood donors assayed in the *Is*-Measles IgG test system and a commercially available test method (see Comparison Testing).

LIMITATIONS

1. The results obtained with the *Is*-Measles IgG Test Kit serve only as an aid to diagnosis and should not be interpreted as diagnostic in themselves.
2. Assay performance characteristics have not been established for visual result determination.
3. The Diamedix *Is*-Measles IgG test system is linear from 20.0 EU/ml (1.0 Index Value) to 100 EU/ml.
4. The test should be performed on serum. The use of whole blood or plasma has not been established.
5. A single positive result only indicates previous immunologic exposure; the level of antibody response or class of antibody may not be used to determine active infection or disease stage.
6. A negative result does not rule out the diagnosis of Rubeola infection. The sample may have been collected before appearance of detectable antibodies. Negative results in suspected early Rubeola infection should be repeated in 4-6 weeks.
7. A significant rise in the level of measles IgG cannot distinguish between primary infection and reinfection with measles. Lack of a significant rise in the level of measles IgG does not exclude the possibility of measles infection.
8. Rare heterotypic responses with rubella virus and varicella virus have been reported for measles virus (7).
9. For individuals experiencing a polyclonal response when infected with a heterotypic virus, a differential diagnosis can be made on the basis of the fact that antibody to the infecting virus type is absent or at a very low titer in the acute-phase specimen, whereas antibody to the viral heterotype is already present (8).
10. The results on serum from immunosuppressed individuals must be interpreted with caution.
11. The performance characteristics of the Diamedix *Is*-Measles IgG Test Kit with automated equipment other than the MAGO[®] Plus Automated EIA Processor, have not been established.

EXPECTED VALUES

The prevalence of measles antibodies can vary depending on age, geographical location, socioeconomic status, race and vaccine usage. Sera from 100 healthy South Florida blood donors (52 female and 48 male) were evaluated in the *Is*-Measles IgG Test Kit. Of the 100 samples, 88 were found to be reactive (positive), 8 were found to be non-reactive (negative) and 4 sera were equivocal. Age distribution, geographic location and prevalence is provided in TABLE 1. Histograms demonstrating the distribution of EU/ml values are shown in FIGURES 1 and 2.

TABLE 1

	Number of donors	Prevalence
Total Number	100	88.0%
Geographic location: South-Eastern US	100	88.0%
Age		
10 – 19	13	84.6%
20 – 29	23	87.0%
30 – 39	40	92.5%
40 – 49	13	84.6%
50 – 59	5	100.0%
60 – 69	6	66.7%

**TABLE 2: Manual
Is-Measles IgG**

		Positive	Negative	Equivocal
Other EIA	Positive	156	10	5
	Negative	1	25	1
	Equivocal	1	5	0
		<i>95% CI*</i>		
Relative Sensitivity		95.6%	91.1 to 98.2	
Relative Specificity		100.0%	86.3 to 100	
Overall Agreement		96.2%	92.3 to 98.5	

FIGURE 1

Is-Measles IgG Reactive Population





FIGURE 2

Is-Measles IgG Non-Reactive Population



PERFORMANCE CHARACTERISTICS

Note: All performance data was generated at Diamedix Corporation, Miami, Florida.

A. Comparison Testing

The Diamedix *Is*-Measles Test Kit results correlated well with results of another commercially available EIA. Sera from 204 normal blood donors were assayed for the presence of Measles IgG antibodies using the Diamedix *Is*-Measles IgG Test Kit and the comparative test. All samples were frozen sera. Sera included 38 samples specifically selected either because they were non-reactive (negative) or had values close to the cut-off by other ELISA methods. Testing was performed both manually and using the MAGO® Plus Automated EIA Processor. Results are shown in Tables 2 and 3.

**TABLE 3: MAGO
Is-Measles IgG**

		Positive	Negative	Equivocal
Other EIA	Positive	152	7	12
	Negative	0	25	2
	Equivocal	1	5	0
		<i>95% CI*</i>		
Relative Sensitivity		94.0%	89.2 to 97.1	
Relative Specificity		96.2%	80.4 to 99.9	
Overall Agreement		94.3%	90.0 to 97.1	

Equivocal results were excluded from calculations.


* Calculated by the Exact Method (9).

NOTE: Please be advised that 'relative' refers to the comparison of the assay's results to that of a similar assay. There was not an attempt to correlate the assay's results with disease presence or absence. No judgment can be made on the comparison assay's accuracy to predict disease.

B. Correlation of Manual and MAGO Plus Results

The Diamedix *Is*-Measles IgG Test Kit has been developed for use both manually and using the MAGO Plus Automated EIA Processor. To further demonstrate the equivalence of the manual and MAGO® procedures, the results of the 204 serum samples tested were plotted and show good correlation (see FIGURE 3).


FIGURE 3



Slope: 0.849
Intercept: 2.18
Corr. Coeff. (r): 0.961
95% CI for r: 0.949 to 0.9708

C. Linearity

Several strongly positive serum specimens were diluted (2-fold) and separate dilutions were assayed in the *Is*-Measles IgG Test Kit both manually and using the MAGO[®] Plus Automated EIA Processor. The titration curve of a representative serum sample is shown in FIGURE 4. The results demonstrate a high degree of linearity throughout the reportable range of the assay.



D. Semi-Quantitative Data

Serum pairs were obtained by preparing multiple serial 2-fold dilutions of several strongly positive sera. Index ratios for dilutions representing a four-fold difference in antibody level were evaluated as a serum pair both manually and using the MAGO[®]. Overall, it was estimated that a 2.1- to a 4.4 -fold (mean 3.2) increase in *Is*-Measles IgG Index value corresponds to a four-fold titer increase in measles IgG antibody level.

E. Cross Reactivity

Twenty sera, non-reactive (negative) for antibodies to measles virus in the *Is*-Measles IgG Test Kit, were tested by ELISA for IgG antibody to varicella, cytomegalovirus, herpes simplex and rubella. 20/20 anti-VZV IgG positive specimens were non-reactive for anti-Rubeola (measles) IgG; 17/17 anti-Rubella IgG positive specimens were non-reactive for anti-Rubeola IgG; 9/9 anti-CMV IgG positive specimens were non-reactive for anti-Rubeola IgG and 12/12 anti-HSV positive specimens were non-reactive for anti-Rubeola IgG. This suggests that no specific cross-reactivity should be expected with the *Is*-Measles IgG Test Kit from these analytes.

F. Precision

Six serum samples (2 negative and 4 positive) plus the *Is*-Measles IgG Test Kit Calibrator, Positive Control and Negative Control, were assayed, in triplicate, in two separate runs on three different days. Testing was performed both manually and using the MAGO[®] Plus Automated EIA Processor. The results for the intra- and interassay precision are shown in Table 4.

TABLE 4

SERUM	Overall MEAN (EU/ml)	MANUAL		MAGO	
		INTRA-CV%	INTER-CV%	INTRA-CV%	INTER-CV%
A (Neg)	7.6	9.9	15.5	5.4	7.0
B (Neg)	4.3	15.1	31.6	10.5	13.2
C (Pos)	32.7	7.5	8.3	6.8	8.0
D (Pos)	50.9	8.6	8.8	5.7	7.5
E (Pos)	85.1	6.3	8.5	8.5	10.2
F (Pos)	103.3	7.2	10.7	5.1	10.6
Calib.	106.4	6.4	8.9	7.2	8.8
Pos. Ctrl.	51.1	6.3	12.1	6.3	9.4
Neg. Ctrl.	3.5	10.5	16.3	16.6	24.3

REFERENCES

1. Snyder, R. C., Gaskins, S. E. and Pieroni, R. E. 1988. Rubeola. *Amer. Fam. Phys.* 137: 175-178.
2. Carter, M. J. and ter Meulen, V. 1987. Measles. In: *Principle and Practice of Clinical Virology*. Zuckerman, A. J., Banatvala, J. E. and Pattison, J. R. (eds). John Wiley and Sons Ltd., New York, p.291-314.
3. Salmi, A. A. Measles Virus. In: *Manual of Clinical Microbiology*. Baron, E.J., Pfaller, M. A., Tenover, F. C. and Tenover, R. H. (eds). 6th Edition, ASM Press, Washington, DC. p. 956-962.
4. CDC Immunization Information - Measles, Centers for Disease Control and Prevention. March 9, 1995.
5. Fuccillo, D. A. and Sever, J. L. 1989. Measles Virus. In: *Diagnostic Procedures for Viral, Rickettsial and Chlamydial infections*. Schmidt, N. J. and Emmons, R. W. (eds). 6th Edition, American Public Health Association, Inc. Washington, DC. p. 713-730.
6. Neumann, P. W., Weber, J. M., Jessamine, A. G. and O'Shaughnessy, M. V. 1985. Comparison of Measles Antihemolysin test, Enzyme-Linked Immunosorbent Assay, and Hemagglutination Inhibition test with Neutralization Test for Determination of Immune Status. *J. Clin. Microbiol.* 22: 296-298.
7. Cremer, N. E. 1984. Anomalous Antibody Responses in Viral Infection: Specific Stimulation or Polyclonal Activation. *J. Clin. Microbiol.* 20: 468-472.
8. Schmidt, N. J. 1985. Varicella Zoster Virus. In: *Manual of Clinical Microbiology*, Lennette, E. H., Balows, A., Hausler, W. J Jr. and Shadomy, H. J. (eds). 4th Edition. American Society for Microbiology, Washington, DC. p.720-727
9. Gardner, M. J. and Altman D. G. 1986. Confidence Intervals rather than P Values: Estimation rather than Hypothesis testing. *Brit. Med. J.* 292: 746-750.
10. Manual Guide – Safety Management No. CDC – 22, “Decontamination of Laboratory Sink Drains to Remove Azide Salts”, Centers for Disease Control and Prevention, Atlanta, GA, April 30, 1976.

Proclin[®] 300 is a registered trademark of Rohm and Haas Corp. Philadelphia, PA.

Diamedix Corporation • A Subsidiary of ERBA Diagnostics, Inc.
14100 NW 57th Court – Miami Lakes, Florida 33014 - USA
(305) 324-2300 / (800) 327-4565
www.erbadiagnostics.com

EC REP Delta Biologicals S.r.l., Via Nicaragua 12/14, 00071 - Pomezia, Rome Italy
Telephone #: +39-06-91190.1 Fax #: +39-069105244



I-720-520
Rev. 8 – June 15

THIS PAGE HAS INTENTIONALLY BEEN LEFT BLANK.