



# Delta Biologicals

A Subsidiary of Erba Diagnostics, Inc.


*Mycobacterium tuberculosis*



*Respiratory Virus*



*Enterovirus*



*Papilloma Virus*




**Microarray for  
multiplex diagnosis**

A MEMBER OF  
**Erba**<sup>®</sup>  
Mannheim

## 8D1082 DR.TBDR/NTM™ IVD Kit


Test for the screening of *Mycobacterium* complex, Rifampicin resistance and identification of 15 non-common tuberculous mycobacterium (NTM) species.



A1, A6, F1, F6	● Hybridization Positive Control
A4	● PCR Positive Control
F4	● Negative control
B1	● <i>M. abscessus</i>
C1	● <i>M. asiaticum</i>
D1	● <i>M. avium</i>
E1	● <i>M. intracellulare</i>
A2	● <i>M. chelonae</i>
B2	● <i>M. fortuitum</i>
C2	● <i>M. goodii</i>
D2	● <i>M. kansasii</i>
E2	● <i>M. lentiflavum</i>
F2	● <i>M. mageritense</i>
A3	● <i>M. marinum</i>
B3	● <i>M. scrofulaceum</i>
C3	● <i>M. shimoidei</i>
D3	● <i>M. szulgii</i>
E3	● <i>M. xenopi</i>
F3	● Mycobacterium
F5	● <i>M. tuberculosis</i> complex
A5	● TB-512WT
B5	● TB-516WT
C5	● TB-523WT
D5	● TB-526WT
E5	● TB-531WT
B6	● TB-mut516a
B4	● TB-mut516b
D4	● TB-mut526a
D6	● TB-mut526b
C4	● TB-mut526c
E4	● TB-mut531a
E6	● TB-mut531b

## 8D1079 DR.HPV-27™ IVD Kit


Test for the genotyping of Papilloma Virus infection (High-risk: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68, 73, 82 - Hypothetical high-risk: 53, 70 - Low-risk: 6, 11, 54, 61, 62, 69, 72, 81, 84)



● Hybridization positive control	● Negative control
● PCR control (β)	● HPV70
● HPV16	● HPV73
● HPV18	● HPV82
● HPV31	● HPV6
● HPV33	● HPV11
● HPV35	● HPV54
● HPV39	● HPV61
● HPV45	● HPV72
● HPV51	● HPV81
● HPV52	● HPV82
● HPV53	● HPV84
● HPV56	● HPV consensus
● HPV58	
● HPV59	
● HPV66	
● HPV68	
● HPV69	

## 8D1040 DR.RV™ IVD (N3) Kit


Test for the detection of 7 major respiratory pathogens in a single assay, including influenza A, influenza B, parainfluenza 1, parainfluenza 2, parainfluenza 3, RSV and adenovirus



A1.A5. E1.E5.	● Hybridization positive control
C3.	● PCR positive control
B2.	● Parainfluenza III
B3.	● Parainfluenza II
B4.	● Parainfluenza I
D1.	● Adenovirus
D2.	● Influenza B
D3.	● Influenza A
D4.	● RSV

## 8D1063 DR.EV™ IVD Kit

Test for the detection of the Enterovirus genus and serotypes enterovirus 71 and coxsackie A16.



A1.A5. E1.E5.	● Hybridization positive control
C3.	● PCR positive control
B2.	● Pan EV
B4.	● EV 71
D2.	● Cox A16



Delta Biologicals Srl - A Subsidiary of Erba Diagnostics, Inc.

Via Nicaragua 12/14

00040 - Pomezia (RM) Italy

Phone: +39 06 91190 Fax: +39 06 9105244

e-mail: info@deltabiologicals.it

F-CHIP rev. 0 10/2012